

Ordering	
CEN 91504	Orbital 310 Centrifuge with Digital Timer and 8 Place Rotor for Open Tube Carriers (Carriers Supplied)
CEN 91505	Orbital 310 Centrifuge with Digital Timer and 6 Place Rotor for Sealed Tube Carriers (Carriers Supplied)

Technical Specifications	
Capacity	8 x 15 ml Open Tube Carriers (CEN 91509) 6 x 15 ml Sealed Tube Carriers (CEN 91510)
Speed	Variable speed control from 1000 to 3200 RPM
RCF (max)	1,700 g
Timer	4 digit, 0-99 mins 59 secs
Analog Displays	Run time remaining Set RPM
Indicators	Min, Sec, Start, Stop 4 Digit LCD time display Run LED
Controls	Digital Timer, Speed Set
End of Cycle	Audible alert
Safety	Lid interlocked power cut off switch
Power	240 V 50 Hz 0.5A
Shipping Weight	16 kg
Unit Dimensions	Height 245 mm, Width 350 mm, Depth 370 mm
Shipping Dimensions	Height 350 mm, Width 400 mm, Depth 450 mm

The Clements range of medical and scientific products is renowned worldwide for quality, reliability and durability

Distributed By:
Bacto Laboratories Pty Ltd 310-312 Elizabeth Dr, Liverpool, 2170 Ph +612 9602 5499 Fax +612 9601 8293 Email: info@bacto.com.au Web Site: www.bacto.com.au

Proudly designed and manufactured in Australia by:
Clements Medical Equipment Pty Ltd 43 Mary Parade (PO Box 137) Rydalmere NSW 2116 Australia Tel +61 2 9638 3285 Fax +61 2 9638 3295 Email: sales@clemequip.com.au Web site: www.clemequip.com.au
 AUSTRALIAN MADE