

Ordering	
CEN 91516	Orbital 325 Centrifuge with Electronic Lid Lock , Digital Timer and 8 Place Rotor for Open Tube Carriers (Carriers Supplied)
CEN 91517	Orbital 325 Centrifuge with Electronic Lid Lock, Digital timer and 6 Place Rotor for Sealed Tube Carriers (Carriers Supplied)

Technical Specifications	
Capacity	8 x 15 ml Open Tube Carriers (CEN 91509) 6 x 15 ml Sealed Tube Carriers (CEN 91510)
Speed	Variable speed control from 1000 to 3200 RPM
RCF (max)	1,700 g
Timer	4 digit 0-99 mins 59 secs
Analog Displays	Set RPM
Indicators	Min, Sec, Start, Stop 4 Digit LCD time display Run LED
Controls	Digital Timer, Speed Set
Safety	Electronic Lid Lock
End of Cycle	Audible alert
Power	240 V 50 Hz 0.5A
Shipping Weight	16 kg
Unit Dimensions	Height 245 mm, Width 350 mm, Depth 370 mm
Shipping Dimensions	Height 350 mm, Width 400 mm, Depth 450 mm

*Specifications subject to change without notice.

The Clements range of medical and scientific products is renowned worldwide for quality, reliability and durability

Distributed By:
Bacto Laboratories Pty Ltd 310-312 Elizabeth Dr, Liverpool, 2170 Ph +612 9602 5499 Fax +612 9601 8293 Email: info@bacto.com.au Web Site: www.bacto.com.au

Proudly designed and manufactured in Australia by:
Clements Medical Equipment Pty Ltd 43 Mary Parade (PO Box 137) Rydalmere NSW 2116 Australia Tel +61 2 9638 3285 Fax +61 2 9638 3295 Email: sales@clemequip.com.au Web site: www.clemequip.com.au
 AUSTRALIAN MADE [®]