

Ordering	
CEN 92115	Orbital 420 Centrifuge with Digital Timer and CEN 91004 Rotor with CEN 91101 (4 x 100 ml) Buckets (Buckets Included)
CEN 92116	Orbital 420 Centrifuge with Digital Timer and CEN 93004 Rotor for Multi Buckets. (Buckets not included)

Technical Specifications	
Speed	Variable speed control from 0 - to 3500 RPM
RCF (max)	2,400 g
Timer	4 digit, 0-99 mins 59 secs
Analog Displays	Run time remaining Set RPM
Indicators	Min, Sec, Start, Stop 4 Digit LCD time display Run LED
End of Cycle	Audible alert
Controls	Digital Timer, Speed Set
Safety	Integrated Electronic Lid Lock
Braking	Manually Operated electric braking
Acceleration	55 to 65 seconds to 3500 rpm
Power	240V 50 Hz 0.5A
Warning Lamps	Imbalance, Power On, Low Brush
Shipping Weight	29 kg
Unit Dimensions	320H • 470W • 500D

*Specifications subject to change without notice.

The Clements range of medical and scientific products is renowned worldwide for quality reliability and durability .

Distributed By:
Bacto Laboratories Pty Ltd 310-312 Elizabeth Dr, Liverpool, 2170 Ph +612 9602 5499 Fax +612 9601 8293 Email: info@bacto.com.au Web Site: www.bacto.com.au

Proudly designed and manufactured in Australia by:
Clements Medical Equipment Pty Ltd 43 Mary Parade (PO Box 137) Rydalmere NSW 2116 Australia Tel +61 2 9638 3285 Fax +61 2 9638 3295 Email: sales@clemequip.com.au Web site: www.clemequip.com.au
 AUSTRALIAN MADE [®]