

The Coombs Test centrifuge is also suitable for blood grouping, typing and cross matching, spinning down samples of whole blood, urines and other body fluids.

Ordering	
CEN 96225	Coombs Test Centrifuge with Digital Timer complete with CEN 96242 Fixed Angle 12 Place Head

Technical Specifications	
Bowl	Stainless steel
Timer	4 digit, 0-99 mins 59 secs
End of Cycle	Audible alert
Indicators	Min, Sec, Start, Stop 4 Digit LCD time display Run LED
Power	220V/240V - 50 Hz.
Speed	2800 rpm
R. C. F.	700 (with Cat. No. CEN 96242 Head)
Safety	Lid interlocked power cut off switch
Unit Dimensions	Height 245 mm, Width 350 mm, Depth 370 mm
Shipping Dimensions	Height 350 mm, Width 400 mm, Depth 450 mm
Shipping Weight	18 kg

*Specifications subject to change without notice.

The Clements range of medical and scientific products is renowned worldwide for quality, reliability and durability

Distributed By:	Proudly designed and manufactured in Australia by:
Bacto Laboratories Pty Ltd 310-312 Elizabeth Dr, Liverpool, 2170 Ph +612 9602 5499 Fax +612 9601 8293 Email: info@bacto.com.au Web Site: www.bacto.com.au	Clements Medical Equipment Pty Ltd 43 Mary Parade (PO Box 137) Rydalmere NSW 2116 Australia Tel +61 2 9638 3285 Fax +61 2 9638 3295 Email: sales@clemequip.com.au Web site: www.clemequip.com.au

