

Ordering	
<b>CEN 96222</b>	Microhaematocrit Centrifuge with 24-place disc head, Digital Timer and Linear Reader.

Technical Specifications	
<b>Capacity</b>	24 x 1.75 mm capillary tubes
<b>Speed</b>	Fixed speed control for 12,500 RPM
<b>RCF (max)</b>	15,500 g
<b>Timer</b>	4 digit, 0 – 99 mins 59 secs
<b>Indicators</b>	Min, Sec, Start, Stop 4 Digit LCD time display Run LED
<b>End of Cycle</b>	Audible alert
<b>Analog Display</b>	Run time remaining
<b>Controls</b>	Digital Timer
<b>Safety</b>	Lid interlocked power cut off switch
<b>Power</b>	240 V 50 Hz 0.5A
<b>Shipping Weight</b>	16 kg
<b>Unit Dimensions</b>	Height 245 mm, Width 350 mm, Depth 370 mm
<b>Shipping Dimensions</b>	Height 350 mm, Width 400 mm, Depth 450 mm

\*Specifications subject to change without notice.

The Clements Linear Reader (CEN 96226) is a compact and lightweight hand held reader. The ergonomic design of the reader ensures operator ease of use and the shaded design of the scale helps to minimise inaccurate readings .

**The Clements range of medical and scientific products is renowned worldwide for quality, reliability and durability**

Distributed By:
Bacto Laboratories Pty Ltd 310-312 Elizabeth Dr, Liverpool, 2170 Ph +612 9602 5499 Fax +612 9601 8293 Email: info@bacto.com.au Web Site: www.bacto.com.au

Proudly designed and manufactured in Australia by:
<b>Clements Medical Equipment Pty Ltd</b> 43 Mary Parade (PO Box 137) Rydalmere NSW 2116 Australia Tel +61 2 9638 3285 Fax +61 2 9638 3295 Email: sales@clemequip.com.au Web site: www.clemequip.com.au
 AUSTRALIAN MADE <sup>®</sup>